


## Basic Heal with Universal viaDNA Remote – Black Box

Select the spooky2 software **Settings** tab



In the **Signal Settings** box type in the **Amplitude** and **Offset** numbers for **Out1** and **Out2**


Signal Settings		Out 1	Out 2	
Duty Cycle		50	50	%
Amplitude		9.9	9.9	V
Offset		-100	-100	%
Phase		0	0	°

In the **Output Control** box ... select **Out2 Follows Program Frequency**

Output Control	
Out 1 =	0 X
Swap Waveform Every	
<input type="radio"/> Out 2 = (Out 1	X
<input type="radio"/> Out 1 fixed at	
<input checked="" type="radio"/> Out 2 Follows Program Frequency	
<input type="radio"/> Out 2 Runs Every Second Frequency	

## Basic Heal with Universal viaDNA Remote – Black Box

In the **Waveform Setup** box ... select **Sine Waveform** for 1 and 2



The oscilloscope box (right hand side of Settings window) should look like this

